


**CONTINUOUS NONINVASIVE  
BLOOD PRESSURE & HEMODYNAMICS**

THE NEXT GENERATION OF  
NONINVASIVE MONITORING



# CONTINUOUS NONINVASIVE HEMODYNAMIC CONTROL

## FULL HEMODYNAMIC PICTURE



### Hemodynamics

- > Continuous noninvasive blood pressure waveform / trendview
- > Cardiac Output  
CO, CI, SV, SI
- > Vascular Resistance  
SVR, SVRI

### Dynamic Fluid Management

- > PPV, SVV

- > Continuous Blood Pressure: Sys, Dia, MAP, Pulse and Upper arm NBP: Sys, Dia

## CONVENIENT CNAP® FINGER SENSOR

-

## NONINVASIVE



### EASY-TO-USE AND QUICK

- > Quick set-up and error-free application
- > Blood pressure waveform and values immediately available

### ACCURATE AND RELIABLE

- > Comparable with invasive clinical standards<sup>1-11, 21</sup>
- > Reliable tracking (e.g. in patients with volatile blood pressure; during Goal Directed Therapy)
- > Noninvasive hemodynamic monitoring can be used as an addition to arterial line

### COST EFFECTIVE

- > Up to 77% cost savings through reusable CNAP® double finger sensor

# EASY-TO-USE QUICK START UP COST EFFECTIVE



## PROVEN ACCURACY IN CLINICAL SETTINGS

- > CNAP® measurements are comparable to invasive arterial line measurements in terms of continuity, accuracy and waveform dynamics.<sup>1,2,3</sup>
- > CNAP® delivers reliable results for the efficient treatment of ICU and ER patients.<sup>4,5, 6,7</sup>
- > CNAP® provides immediate hemodynamic status and detects blood pressure drops during the induction of anesthesia.<sup>8</sup>
- > CNAP® shows outstanding performance in monitoring pediatric patients without an arterial catheter.<sup>9,10</sup>

## EASY & RELIABLE TOOL FOR RESEARCH

- > Noninvasive measurement
- > Easy-to-use: all from *one* sensor
- > Reliability clinically validated

## FAST & ACCURATE HEMODYNAMIC OVERVIEW

- > Early recognition<sup>15</sup>
- > Fast intervention
- > Detection of hemodynamic reactions
- > ...without arterial catheter

## REDUCING RISK & IMPROVING OUTCOME THROUGH GOAL DIRECTED THERAPY

- > Noninvasive CNAP® PPV / SVV is an accurate predictor of fluid responsiveness in anaesthetized patients.<sup>18,19</sup>
- > Goal directed therapy with CNAP® HD significantly reduces postoperative infections, organ complications and number of transfusions.<sup>20</sup>
- > Noninvasive CO with CNAP® HD performs comparably to invasive CO monitoring.<sup>21</sup>

**“Given the fact that CNAP® is a reliable device to assess the arterial AP continuously, [...] its noninvasiveness facilitates its use for any operation with a need to assess, document, and maintain hemodynamic stability.”<sup>1</sup>**

**“CNAP® can be used as an alternative to intra-arterial pressure”<sup>4</sup>**



- 1 Jeleazcov, C. et al. Precision and accuracy of a new device (CNAP®) for continuous noninvasive arterial blood pressure monitoring: assessment during general anaesthesia. *BJA*. 105(3):264-272 (2010).
- 2 Ilies, C., Investigation of the agreement of a continuous non-invasive arterial pressure device in comparison with invasive radial artery measurement. *BJA*. 108(2):202-10. doi: 10.1093/bja/aer394 (2012).
- 3 Blais, M. et al. Continuous non-invasive arterial pressure measurement: Evaluation of CNAP™ device during vascular surgery. *Ann Fr Anesth Reanim*, doi:10.1016/j.annfar.2010.05.002 (2010)
- 4 Jagadeesh, AM., A comparison of a continuous noninvasive arterial pressure (CNAP™) monitor with an invasive arterial blood pressure monitor in the cardiac surgical ICU. *Ann Card Anaesth*. Jul-Sep;15(3):180-4. doi: 10.4103/09719784.97973 (2012).
- 5 Ilies, C. et al. Comparison of a continuous noninvasive arterial pressure device with invasive measurements in cardiovascular postsurgical intensive care patients: A prospective observational study. *European Journal of Anaesthesiology*, 31, 1–9. doi:10.1097/EJA.0000000000001366 (2014).
- 6 Wagner, J. Y. et al. Noninvasive continuous versus intermittent arterial pressure monitoring: evaluation of the vascular unloading technique (CNAP device) in the emergency department. *Scandinavian Journal of Trauma, Resuscitation and Emergency Medicine*, 22(1), 8. doi:10.1186/1757-7241-22-8 (2014).
- 7 Wagner, J. Y. et al. Continuous noninvasive arterial pressure measurement using the volume clamp method : an evaluation of the CNAP device in intensive care unit patients. *J Clin Monit Comput*, online. doi:10.1007/s10877-015-9670-2 (2015).
- 8 Kumar, A., et al. Evaluation of continuous non - invasive arterial pressure monitoring during induction of general anaesthesia in patients undergoing cardiac surgery. *Indian J Anaesth*, 59(1), 21–25. doi:10.4103/0019-5049.149444 (2015).
- 9 Kako, H. et al. Accuracy of the CNAP™ monitor, a noninvasive continuous blood pressure device, in providing beat-to-beat blood pressure readings in pediatric patients weighing 20–40 kilograms. *Paediatric Anaesthesia*, 1–5. doi:10.1111/pan.12173 (2013).
- 10 Dewhirst, E. et al. Accuracy of the CNAP monitor, a noninvasive continuous blood pressure device, in providing beat-to-beat blood pressure readings in the prone position. *Journal of Clinical Anesthesia*, 1–4. doi:10.1016/j.jclineane.2013.01.01(2013).
- 11 Gonzales, J. U. et al. Arterial stiffness is higher in older adults with increased perceived fatigue and fatigability during walking. *Experimental Gerontology*. doi:10.1016/j.exger.2014.12.005 (2014).
- 12 Lee JF, et al. The magnitude of heat-stress induced reductions in cerebral perfusion does not predict heat-stress induced reductions in tolerance to a simulated hemorrhage. *Journal of Applied Physiology*, 114(1), 37–44. (2013).
- 13 Sng, B. L. et al. Closed-loop double-vasopressor automated system vs manual bolus vasopressor to treat hypotension during spinal anaesthesia for caesarean section: a randomised controlled trial. *Anaesthesia*, 1–9. doi:10.1111/anae.12460 (2013).
- 14 Cornick, J. E. et al. Consequences of objective self-awareness during exercise. *Health Psychology Open*, 2(2), 2055102915598088. doi:10.1177/2055102915598088 (2015).
- 15 Benes, J., et al. Continuous non-invasive monitoring improves blood pressure stability in upright position: randomized controlled trial. *Journal of Clinical Monitoring and Computing*. doi:10.1007/s10877-014-9586-2 (2014).
- 16 Ilies, C. et al. Detection of hypotension during Caesarean section with continuous non-invasive arterial pressure device or intermittent oscillometric arterial pressure measurement. *British Journal of Anaesthesia*, 3–9. doi:10.1093/bja/aes224 (2012).
- 17 Siebig, S. et al. Continuous non-invasive arterial pressure technique improves patient monitoring during interventional endoscopy. *International Journal of Medical Sciences*, 6(1), 37–42. Retrieved from <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=2631161&tool=pmcentrez&rendertype=abstract> (2009).
- 18 Blais, M. et al., The ability of pulse pressure variations obtained with CNAP™ device to predict fluid responsiveness in the operating room; *Anesthesia and analgesia*, 523-28 (2011).
- 19 Monnet, X. et al., Prediction of fluid responsiveness by a continuous non-invasive assessment of arterial pressure in critically ill patients: comparison with four other dynamic indices. *British Journal of Anaesthesia* (2012).
- 20 Benes, J et al. Fluid management guided by a CNAP device is associated with decreased postoperative morbidity after total knee and hip replacement. *BMC Anesthesiology*, 15(1), 148 (2015).
- 21 Wagner, J. Y., et al. A comparison of volume clamp method-based continuous noninvasive cardiac output (CNCO) measurement versus intermittent pulmonary artery thermodilution in postoperative cardiothoracic surgery patients. *Journal of Clinical Monitoring and Computing*, 1–10. <http://doi.org/10.1007/s10877-017-0027-x> (2017).

## TECHNICAL SPECIFICATIONS

CNAP® – CONTINUOUS NONINVASIVE ARTERIAL PRESSURE			NBP – OSCILLOMETRIC BLOOD PRESSURE	
Measuring range	Sys:	40 - 250 mmHg	Measuring range	Sys: Adult 40 - 260 mmHg, Pediatric 40 - 230 mmHg
	Dia:	30 - 210 mmHg		Dia: Adult 20 - 200 mmHg, Pediatric 20 - 160 mmHg
	Mean:	35 - 230 mmHg		
	Pulse rate:	30 - 200 bpm		
Degree of protection	BF (defibrillation proof)		Degree of protection	BF (defibrillation proof)
Automatic scaling to brachial pressure (NBP)				
CNAP® HEMODYNAMICS: CO, CI, SV, SVR, SVI, SVRI				
Measuring range	CO	0.0 - 15 l/min	CI	0.0 - 8 l/min/m <sup>2</sup>
	SV	0 - 200 ml	SI	0 - 100 ml/m <sup>2</sup>
	SVR	0 - 3000 dyne*s/cm <sup>5</sup>	SVRI	0 - 6000 dyne*s*m <sup>2</sup> /cm <sup>5</sup>
FLUID RESPONSIVENESS: CNAP® PPV AND SVV				
Measuring range	PPV: 0.2 - 40%;		SVV: 0 - 40%	
ELECTRICAL				
Nominal voltage	100 - 240 VAC		Battery:	sealed lead-gel, operating time: 2 hours (fully charged battery)
Supply frequency	~50/60 Hz			
PHYSICAL				
Weight	7,5 kg (16,6 lbs) including accessories and cables			
Height	280 x 270 x 250 mm (11 x 10,6 x 9,8 inch)			
ENVIRONMENTAL				
Temperature	operation:	10°C - 40°C (50°F - 104°F)	storage:	0°C - 40°C (32°F - 104°F)
Humidity	operation:	15% - 85% non condensing	storage:	15% - 95%, non condensing, wrapped
Altitude	operation:	647 - 1060 hPa	storage:	500 - 1060 hPa
SCREEN				
Type	TFT-LCD, 800 x 600 pixel			
Size	8,4 inch diagonally			
USER INTERFACE				
Controls	click-wheel control, fast access keys			
Indicators	visual and audible alarm indication, battery status, printer status, power LED			
Trend Display	customized configuration: numeric, graphic, alarm history			
ADJUSTABLE ALARMING SYSTEM				
Alarms	physiological: med priority, technical: low priority			
CONNECTIVITY				
BP Wave Out	easy integration in all standard patient monitoring systems (2 - 10 VDC supply voltage)			
AUX Analog Out	analog output of calibrated continuous blood pressure waveform and additional configurable hemodynamic parameters (0V to 5V; -5V to 5V)			
USB PORT				
Version	USB 1.1 (bandwidth 12 Mbits/s)			
PRINTER				
Type	integrated thermal printer, 58 mm			
COMPLIANCE AND APPROVALS				
Safety class II (IEC 60601)	> IEC 60601-1	> IEC 60601-1-6	> EN 1060-4 (NBP)	
Class II b (93/42/EEC)	> IEC 60601-1-2	> IEC 60601-1-8	> ISO 81060-2 (NBP)	
Patient applied part type BF	(defibrillation proof)	> IEC 80601-2-30		
INTELLECTUAL PROPERTY				
Patents	> US 6,669,648	> US 8,343,062	> JP 20075508872	
	> EP 1 179 991	> EP 2 493 370	> CN 102647940	
	> US 8,114,025	> US 8,814,800 B2		
	> EP 1 675 507	> EP 2 493 373	plus additional 66 patents	

The CNAP® Monitor is CE approved. All parameters in section "CNAP® hemodynamics" and "fluid responsiveness" currently have no FDA clearance.

## CNAP® – Setting new standards for continuous and noninvasive hemodynamic monitoring.



local distributor:

CNSystems Medizintechnik GmbH  
Reininghausstrasse 13, 8020 Graz/Austria  
Tel.: +43 316 723456-0, Fax: +43 316 723456-2  
E-Mail: office@cnsystems.com, www.cnsystems.com