


Extension Cables










The patient extension cables are designed to connect to most external cardiac pacemakers currently available. Each cable is 2.5 metres long, autoclavable and supplied with either a 2mm collet terminal or a Universal Connector. Standard 2mm pins and Heart Wires can be attached to the collet terminals, whilst the Universal Connector will also accept shrouded 2mm pins.

The threshold cables are designed for use with a pacing analyser. All cables are 3 metres long and have colour-coded, insulated alligator clips for connection to lead connector pins, clamps, wires or patient tissue. The 4220 is supplied sterile for single use only. All other models are autoclavable.

Atrium (A) or Ventricular (V) indicators may be specified on selected models.

Extension Cables


Model Reference Chart

| | Proximal Connection | | | | Distal Connection | | |
|-----------------------|---|---|---|---|--|---|---|
| |  |  |  |  |  |  |  |
| Patient Cables | | | | | | | |
| E4231 | ✓ | | | | ✓ | | |
| E4231U | ✓ | | | | | ✓ | |
| E4233 | | ✓ | | | ✓ | | |
| E4233U | | ✓ | | | | ✓ | |
| E4235* | | | ✓ | | ✓ | | |
| 4265* | | | ✓ | | | ✓ | |
| 5235 | | | | ✓ | ✓ | | |
| 5265 | | | | ✓ | | ✓ | |

Threshold Cables

| | | | | | | | |
|--------|---|--|---|--|--|--|---|
| 4205 | ✓ | | | | | | ✓ |
| 4280* | | | ✓ | | | | ✓ |
| 4220** | ✓ | | | | | | ✓ |

Adaptor
4810



Adapts standard 2mm collet terminals with fixed 1.25" spacing to accept catheters supplied with 2mm shrouded pins

Other lengths and proximal / distal connector combinations available on request.

* Available with Atrium (A) or Ventricular (V) indicators at each end of the cable

** Supplied sterile for Single use only

APC Cardiovascular Ltd.
18 Macon Court, Herald Drive,
Crewe, Cheshire CW1 6EA, England
Tel: +44 (0)1270 216142
Fax: +44 (0)1270 216251
Email: info@apc-cardiovascular.co.uk

APC Cardiovascular Ltd.
International Branch Office Nederland
Blauwwater 23-a, 5951 DB, Belfeld, Nederland
Tel: +31 (0)77 475.48.91
Fax: +31 (0)77 475.48.93
Email: apcni@apc-cardiovascular.co.uk

www.apc-cardiovascular.co.uk