


# Myo-Pace™

## Unipolar Myocardial Heart Wires



The Myo-Pace range of unipolar myocardial heart wires is designed specifically for use as a temporary cardiac pacing wire after open heart surgery.

The polyethylene insulation provides the smooth surface necessary to minimise tissue trauma both on placing and withdrawal of the heart wire, whilst the multifilament braided stainless steel wire construction ensures maximum flexibility and strength. Break-off needles are available to permit easy connection to the external pacemaker.

# Myo-Pace™

## Unipolar Myocardial Heart Wires

- Braided stainless steel wire for strength and flexibility
- Smooth polyethylene insulation to minimise tissue trauma
- Spread anchor models available to provide more stable placement
- Colour coded insulation enables easy identification of implanted wires
- Break-off needles permit easy connection to the external pacemaker

Model Number	Colour	Size		Length	Needle*
		USP	Metric		
<b>Standard Needle</b>					
3140	Blue	3-0	2	60cm	HR 26 + GS 63
3145	Transparent	3-0	2	60cm	HR 26 + GS 63
3150	Blue	3-0	2	60cm	HR 17 + GS 63
3155	Transparent	0	3.5	60cm	HR 26 + GS 63

### Break Off Needle

3140B	Blue	3-0	2	60cm	HR 26 + GS 63
3145B	Transparent	3-0	2	60cm	HR 26 + GS 63
3155B	Transparent	0	3.5	60cm	HR 26 + GS 63

### Break Off Needle With Spread Anchor

3170	Blue	3-0	2	60cm	HR 26 + GS 63
3175	Transparent	3-0	2	60cm	HR 26 + GS 63

### Specifications

Electrode:	Multifilament stainless steel wire
Insulation:	Polyethylene
Needles:	Stainless Steel
Sterilisation Process:	Ethylene Oxide (Eto)

### Ordering Information

Myo-Pace myocardial heart wires are packaged 36 units per box

Other sizes and lengths are available on request.

\*HR 17 = half circle round bodied needle, length: 17mm

\*HR 26 = half circle round bodied needle, length: 26mm

\*GS 63 = straight cutting needle, length: 63mm

APC Cardiovascular Ltd.  
18 Macon Court, Herald Drive,  
Crewe, Cheshire CW1 6EA, England  
Tel: +44 (0)1270 216142  
Fax: +44 (0)1270 216251  
Email: info@apc-cardiovascular.co.uk

APC Cardiovascular Ltd.  
International Branch Office Nederland  
Blauwwater 23-a, 5951 DB, Belfeld, Nederland  
Tel: +31 (0)77 475.48.91  
Fax: +31 (0)77 475.48.93  
Email: apcni@apc-cardiovascular.co.uk

[www.apc-cardiovascular.co.uk](http://www.apc-cardiovascular.co.uk)