

Cable Overview

Cable connections for external devices

Cardiac Rhythm Management

External Devices

Cable Overview

Cable Overview

Cable Overview

Cable connections for external devices



© BIOTRONIK SE & Co. KG
All rights reserved. Specifications
are subject to modification,
revision and improvement.

BIOTRONIK SE & Co. KG
Voermannstraße 1
12359 Berlin Germany
Tel. +49 (0)30 68905-0
Fax +49 (0)30 6852804
sales@biotronic.com
www.biotronic.com

 **BIOTRONIK**
excellence for life

 **BIOTRONIK**

 **BIOTRONIK**
excellence for life

1. Reocor D

2. Reocor S

3. Reliaty

4. Renamic PSA/Renamic

5. EDP 20/A

6. EDP 30/A + EDP 30/S








7. EDP 20/B

8. EDP 30/B + EDP 30/BP

9. ERA 300/ERA 3000











10. ICS 3000 IM/ICS 3000

Temporary catheter with shrouded 2 mm pins

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2×						
 2×	 TC-Adapt 4×	 PK-175			- PK-175: 333 959 - TC-Adapt: 121 340 (set of 2 adapters)	1, 2
 2×	 TC-Adapt 4×	 PK-83 2×			- PK-83 (1.5 m): 128 562 - PK-83 (2.5 m): 128 563 - TC-Adapt: 121 340 (set of 2 adapters)	1 1
 2×	 PA-1-B	 PK-67			- PK-67-S (0.8 m): 128 085 - PK-67-L (2.6 m): 123 672 - PA-1-B: 123 751	1 1 1

- 121°C; 1.1 bar; 20 min
- 134°C; 2.2 bar; 5 min
- Single-use only

Reocor D External pacemaker (dual-chamber)

Heartwire (temporary pacing lead) with shrouded 2 mm adapter*						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
						
	 TC-Adapt 4x	 PK-175			• PK-175: 333 959 • TC-Adapt: 121 340 (set of 2 adapters)	1, 2
	 TC-Adapt 4x	 PK-83 2x			• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1 1
	 PA-1-B	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-B: 123 751	1 1 1















1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min
3 Single-use only

* This table shows cable connections for dual chamber devices, cable connections for single chamber devices see Reocor S

Reocor D External pacemaker (dual-chamber)

1

Heartwire (temporary pacing lead) with break-off needle or flexible tip (max. diameter: 2.3 mm)*

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2x		 PK-175			• PK-175: 333 959	1, 2
 2x		 PK-83 2x			• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563	1 1
 2x		 PK-141			• PK-141: 353 182	1, 2
 2x	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1 1 1
 2x	 PK-155 2x	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1 1 3
 2x		 PK-82 2x			• PK-82: 128 564	1









1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min

3 Single-use only

* This table shows cable connections for dual chamber devices, cable connections for single chamber devices see Reocor S

Reocor D External pacemaker (dual-chamber)










Implantable leads*						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-82 2x			• PK-82: 128 564	1
		 PK-141			• PK-141: 353 182	1, 2
	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
					• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
	 PK-155 2x	 PK-67			• PK-155: 337 358	3

- 1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min
3 Single-use only

* This table shows cable connections for dual chamber devices, cable connections for single chamber devices see Reocor S

Reocor S External pacemaker (single-chamber)

Temporary catheter with shrouded 2 mm pins














Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
						
	 TC-Adapt 2x	 PK-83			• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1
	 TC-Adapt 2x	 PK-83-B			• PK-83-B (1.5 m): 347 606 • PK-83-B (2.5 m): 347 485 • TC-Adapt: 121 340 (set of 2 adapters)	1
	 PA-1-C	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-C: 349 723	1
	 PA-1-B*	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-B: 123 751	1

- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

Reacor S External pacemaker (single-chamber)

Heartwire (temporary pacing lead) with shrouded 2 mm adapter

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2x						
 2x	 TC-Adapt 2x	 PK-83			• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1
 2x	 TC-Adapt 2x	 PK-83-B			• PK-83-B (1.5 m): 347 606 • PK-83-B (2.5 m): 347 485 • TC-Adapt: 121 340 (set of 2 adapters)	1
 2x	 PA-1-C	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-C: 349 723	1
 2x	 PA-1-B*	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-B: 123 751	1













1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min

3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

Heartwire (temporary pacing lead) with break-off needle or flexible tip (max. diameter: 2.3 mm)

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
		 PK-83			• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563	1 1
		 PK-83-B			• PK-83-B (1.5 m): 347 606 • PK-83-B (2.5 m): 347 485	1 1
		 PK-141*			• PK-141: 353 182	1, 2
	 PA-4*	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1 1 1
	 PK-155*	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1 1 3
		 PK-82			• PK-82: 128 564	1
		 PK-175*			• PK-175: 333 959	1, 2







1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min

3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

Reocor S External pacemaker (single-chamber)

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2*	 PK-67			• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
					• PA-2: 123 157	1
		 PK-82			• PK-82: 128 564	1
		 PK-141*			• PK-141: 353 182	1, 2
	 PA-4*	 PK-67			• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
					• PA-4: 123 090	1
	 PK-155	 PK-67*			• PK-67-S (0.8 m): 128 085	1, 3
					• PK-67-L (2.6 m): 123 672	1
					• PK-155: 337 358	3








1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min








3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

Implantable leads

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141		 A RV LV	• PK-141: 353 182	1, 2
	 PA-4	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1, 3
	 PA-2	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141		 A RV LV	• PK-141: 353 182	1, 2
	 PA-4	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155 2x	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1, 3

- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085	1
	 PA-2*	 PK-67			• PK-67-L (2.6 m): 123 672	1
	 PK-141	 PK-141		 A RV LV	• PK-141: 353 182	1, 2
		 PK-141*				
	 PA-4	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085	1
	 PA-4*	 PK-67			• PK-67-L (2.6 m): 123 672	1
	 PK-155 3x	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085	1
		 PK-67*			• PK-67-L (2.6 m): 123 672	1
					• PK-155: 337 358	3










1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min

3 Single-use only

* Using a dual chamber cable for a single/three chamber device, do always connect at the ventricular cable pin

Renamic PSA Pacing system analyzer (programmer module)

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141			• PK-141: 353 182	1, 2
	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1
	 PA-2	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141			• PK-141: 353 182	1, 2
	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155 2x	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	3









1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min


3 Single-use only

* Using a dual chamber cable for a single/three chamber device, do always connect at the ventricular cable pin

Renamic PSA Pacing system analyzer (programmer module)

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085	1
	 PA-2*				• PK-67-L (2.6 m): 123 672	1
	 PA-4	 PK-141		 A RV LV	• PK-141: 353 182	1, 2
		 PK-141*				
	 PA-4*	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085	1
	 PA-4*				• PK-67-L (2.6 m): 123 672	1
	 PK-155 3x	 PK-67		 A RV LV	• PK-67-S (0.8 m): 128 085	1
		 PK-67*			• PK-67-L (2.6 m): 123 672	1
					• PK-155: 337 358	3










Renamic Programmer

ECG cable						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 Electrode Clip (set of 4 units)	 PK-222 4x			• PK-222 (3 m): 335 284 • PK-Electrode Clip: 340 293 (set of 4 units)	No











- 121°C; 1.1 bar; 20 min
- 134°C; 2.2 bar; 5 min
- Single-use only

* Using a dual chamber cable for a single/three chamber device, do always connect at the ventricular cable pin

Temporary catheter with shrouded 2 mm pins






Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 TC-Adapt 2x				• TC-Adapt: 121 340 (set of 2 adapters)	
	 TC-Adapt 2x	 PK-83	 TC-Adapt 2x		• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1 1
	 TC-Adapt 2x	 PK-82	 TC-Adapt 2x		• PK-82: 128 564 • TC-Adapt: 121 340 (set of 2 adapters)	1
	 TC-Adapt 2x	 PK-81	 TC-Adapt 2x		• PK-81 (1.5 m): 128 566 • TC-Adapt: 121 340 (set of 2 adapters)	1

Heartwire (temporary pacing lead) with 2 mm adapter



Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2x	 TC-Adapt 2x				• TC-Adapt: 121 340 (set of 2 adapters)	
 2x	 TC-Adapt 2x	 PK-83	 TC-Adapt 2x		• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1 1
 2x	 TC-Adapt 2x	 PK-81	 TC-Adapt 2x		• PK-81 (1.5 m): 128 566 • TC-Adapt: 121 340 (set of 2 adapters)	1

EDP 20/A External pacemaker (single-chamber)

Heartwire (temporary pacing lead) with break-off needle or flexible tip (max. diameter: 2.3 mm)















Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
		 PK-83	 TC-Adapt 2x		• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1
		 PK-82	 TC-Adapt 2x		• PK-82: 128 564 • TC-Adapt: 121 340 (set of 2 adapters)	1

Implantable leads

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
		 PK-82	 TC-Adapt 2x		• PK-82: 128 564 • TC-Adapt: 121 340 (set of 2 adapters)	1

EDP 30/A + EDP 30/S External pacemaker (dual-chamber)

Temporary catheter with shrouded 2 mm pins

















Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2 ×	 TC-Adapt 4 ×				TC-Adapt: 121 340 (set of 2 adapters)	
 2 ×	 TC-Adapt 4 ×	 PK-83 2 ×	 TC-Adapt 4 ×		- PK-83 (1.5 m): 128 562 - PK-83 (2.5 m): 128 563 - TC-Adapt: 121 340 (set of 2 adapters)	1 1
 2 ×	 TC-Adapt 4 ×	 PK-175 + VK-99			- PK-175: 333 959 - VK-99: 333 960 - TC-Adapt: 121 340 (set of 2 adapters)	1, 2 3
 2 ×	 TC-Adapt 4 ×	 PK-81 2 ×	 TC-Adapt 4 ×		- PK-81 (1.5 m): 128 566 - TC-Adapt: 121 340 (set of 2 adapters)	1
 2 ×	 TC-Adapt 4 ×	 PK-82 2 ×	 TC-Adapt 4 ×		- PK-82: 128 564 - TC-Adapt: 121 340 (set of 2 adapters)	1

1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min









EDP 30/A + EDP 30/S External pacemaker (dual-chamber)

Heartwire (temporary pacing lead) with shrouded 2 mm adapter

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 TC-Adapt 4x				• TC-Adapt: 121 340 (set of 2 adapters)	
	 TC-Adapt 4x	 PK-83 2x	 TC-Adapt 4x		• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1 1
	 TC-Adapt 4x	 PK-175 + VK-99			• PK-175: 333 959 • VK-99: 333 960 • TC-Adapt: 121 340 (set of 2 adapters)	1, 2 3
	 TC-Adapt 4x	 PK-81 2x	 TC-Adapt 4x		• PK-81 (1.5 m): 128 566 • TC-Adapt: 121 340 (set of 2 adapters)	1
	 TC-Adapt 4x	 PK-82 2x	 TC-Adapt 4x		• PK-82: 128 564 • TC-Adapt: 121 340 (set of 2 adapters)	1




1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min

EDP 30/A + EDP 30/S External pacemaker (dual-chamber)

Heartwire (temporary pacing lead) with break-off needle or flexible tip (max. diameter: 2.3 mm)						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2×		 PK-83 2×	 TC-Adapt 4×		• PK-83 (1.5 m): 128 562 • PK-83 (2.5 m): 128 563 • TC-Adapt: 121 340 (set of 2 adapters)	1 1
 2×		 PK-175 + VK-99			• PK-175: 333 959 • VK-99: 333 960	1, 2 3
 2×		 PK-82 2×	 TC-Adapt 4×		• PK-82: 128 564 • TC-Adapt: 121 340 (set of 2 adapters)	1

1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min
3 Single-use only











EDP 30/A + EDP 30/S External pacemaker (dual-chamber)

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
		 PK-82 2x	 TC-Adapt 4x		- PK-82: 128 564 - TC-Adapt: 121 340 (set of 2 adapters)	1

1 121°C; 1.1 bar; 20 min

EDP 20/B External pacemaker (single-chamber)

Temporary catheter with shrouded 2 mm pins

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 TC-Adapt 2x	 PK-83-B			• PK-83-B (1.5 m): 347 606 • PK-83-B (2.5 m): 347 485 • TC-Adapt: 121 340 (set of 2 adapters)	1
	 PA-1-C	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-C: 349 723	1
	 PA-1-B*	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-B: 123 751	1
	 TC-Adapt 2x	 PK-141*			• PK-141: 353 182 • TC-Adapt: 121 340 (set of 2 adapters)	1, 2
	 TC-Adapt + PA-4* 2x	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090 • TC-Adapt: 121 340 (set of 2 adapters)	1
	 TC-Adapt + PK-155 2x	 PK-67*			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358 • TC-Adapt: 121 340 (set of 2 adapters)	3

1 121°C; 1.1 bar; 20 min

















2 134°C; 2.2 bar; 5 min

3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

EDP 20/B External pacemaker (single-chamber)

Heartwire (temporary pacing lead) with shrouded 2 mm adapter

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2x	 TC-Adapt 2x	 PK-83-B			- PK-83-B (1.5 m): 347 606 - PK-83-B (2.5 m): 347 485 - TC-Adapt: 121 340 (set of 2 adapters)	1
 2x	 PA-1-C	 PK-67			- PK-67-S (0.8 m): 128 085 - PK-67-L (2.6 m): 123 672 - PA-1-C: 349 723	1
 2x	 PA-1-B*	 PK-67			- PK-67-S (0.8 m): 128 085 - PK-67-L (2.6 m): 123 672 - PA-1-B: 123 751	1
 2x	 TC-Adapt 2x	 PK-141*			- PK-141: 353 182 - TC-Adapt: 121 340 (set of 2 adapters)	1, 2
 2x	 TC-Adapt + PA-4* 2x	 PK-67			- PK-67-S (0.8 m): 128 085 - PK-67-L (2.6 m): 123 672 - PA-4: 123 090 - TC-Adapt: 121 340 (set of 2 adapters)	1
 2x	 TC-Adapt + PK-155 2x	 PK-67*			- PK-67-S (0.8 m): 128 085 - PK-67-L (2.6 m): 123 672 - PK-155: 337 358 - TC-Adapt: 121 340 (set of 2 adapters)	3

1 121°C; 1.1 bar; 20 min





2 134°C; 2.2 bar; 5 min

3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

EDP 20/B External pacemaker (single-chamber)






Heartwire (temporary pacing lead) with break-off needle or flexible tip (max. diameter: 2.3 mm)

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
		 PK-83-B			• PK-83-B (1.5 m): 347 606 • PK-83-B (2.5 m): 347 485	1 1
		 PK-141*			• PK-141: 353 182	1, 2
	 PA-4*	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1 1 1
	 PK-155	 PK-67*			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1 1 3

- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin

EDP 20/B External pacemaker (single-chamber)














Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2*	 PK-67			• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
	 PA-4*	 PK-141*			• PK-141: 353 182	1, 2
	 PK-155	 PK-67			• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
	 PK-155	 PK-67*			• PA-4: 123 090	1
					• PK-67-S (0.8 m): 128 085	1
	 PK-155	 PK-67*			• PK-67-L (2.6 m): 123 672	1
					• PK-155: 337 358	3

1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min
3 Single-use only

* Using a dual chamber cable/adapter for a single chamber device, do always connect at the ventricular cable pin





EDP 30/B + EDP 30/BP External pacemaker (dual-chamber)

Temporary catheter with 2 mm pins

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2×	 TC-Adapt 4×	 PK-175			• PK-175: 333 959 • TC-Adapt: 121 340 (set of 2 adapters)	2
 2×	 PA-1-B	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-B: 123 751	1 1 1
 2×	 TC-Adapt 4×	 PK-141			• PK-141: 353 182 • TC-Adapt: 121 340 (set of 2 adapters)	1, 2
 2×	 TC-Adapt + PA-4 4×	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090 • TC-Adapt: 121 340 (set of 2 adapters)	1 1 1
 2×	 TC-Adapt + PK-155 4× 2×	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358 • TC-Adapt: 121 340 (set of 2 adapters)	1 1 3

- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only









EDP 30/B + EDP 30/BP External pacemaker (dual-chamber)

Heartwire (temporary pacing lead) with shrouded 2 mm adapter						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 4 x	 TC-Adapt 4 x	 PK-175			• PK-175: 333 959 • TC-Adapt: 121340 (set of 2 adapters)	1, 2
 4 x	 PA-1-B	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-1-B: 123751	1 1 1

1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min
3 Single-use only


EDP 30/B + EDP 30/BP External pacemaker (dual-chamber)

Heartwire (temporary pacing lead) with break-off needle or flexible tip (max. diameter: 2.3 mm)

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
 2×		 PK-175			• PK-175: 333 959	1, 2
 2×		 PK-141			• PK-141: 353 182	1, 2
 2×	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1 1 1
 2×	 PK-155 2×	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1 1 3







- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only

EDP 30/B + EDP 30/BP External pacemaker (dual-chamber)





Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67			• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
	 PA-2	 PK-141			• PK-141: 353 182	1, 2
	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085	1
					• PK-67-L (2.6 m): 123 672	1
	 PK-155 2x	 PK-67			• PA-4: 123 090	1
					• PK-67-S (0.8 m): 128 085	1
	 PK-155 2x	 PK-67			• PK-67-L (2.6 m): 123 672	1
					• PK-155: 337 358	3

1 121°C; 1.1 bar; 20 min
2 134°C; 2.2 bar; 5 min
3 Single-use only

Implantable leads

Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141			• PK-141: 353 182	1, 2
	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1, 3
	 PA-2	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141			• PK-141: 353 182	1, 2
	 PA-4	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155 2x	 PK-67			• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1 1 3










- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141		 RV LV	• PK-141: 353 182	1, 2
	 PA-4	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1
	 PA-2	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-2: 123 157	1
		 PK-141		 RV LV	• PK-141: 353 182	1, 2
	 PA-4	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PA-4: 123 090	1
	 PK-155 2x	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085 • PK-67-L (2.6 m): 123 672 • PK-155: 337 358	1
	 PA-3	 PK-144		 RV LV	• PK-144: 353 906 • PA-3: 353 969	1


- 1 121°C; 1.1 bar; 20 min
- 2 134°C; 2.2 bar; 5 min
- 3 Single-use only

* Using a dual chamber cable for a single/three chamber device, do always connect at the ventricular cable pin

ICS 3000 IM Programmer with implant module

Implantable leads						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 PA-2	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085	1
	 PA-2*				• PK-67-L (2.6 m): 123 672	1
	 PA-2*	 PK-141		 RV LV	• PK-141: 353 182	1, 2
		 PK-141*				
	 PA-4	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085	1
	 PA-4*				• PK-67-L (2.6 m): 123 672	1
	 PA-4*	 PK-67		 RV LV	• PK-67-S (0.8 m): 128 085	1
		 PK-67*			• PK-67-L (2.6 m): 123 672	1
	 PK-155 3x	 PK-67		 RV LV	• PK-155: 337 358	3

ICS 3000 / ICS 3000 IM

ECG cable						
Electrode	Patient side	BIOTRONIK cable	Adapter	Device side	Order number	Sterilization
	 Electrode Clip (set of 4 units)	 PK-222 4x			• PK-222 (3 m): 335 284 • PK-Electrode Clip: 340 293 (set of 4 units)	3

1 121°C; 1.1 bar; 20 min

2 134°C; 2.2 bar; 5 min

3 Single-use only

* Using a dual chamber cable for a single/three chamber device, do always connect at the ventricular cable pin