

TEMPORARY BIPOLAR PACING CATHETER

TEMPORARY BIPOLAR PACING CATHETER STANDARD OR SOFT TIP – TORQUE DIRECTED

- intravenous electrode for pre- and postoperative and emergency stimulation
- interim support during pacemaker exchange
- intracardial monitoring
- soft tip version with flexible distal end for atraumatic placement in delicate vessel structures

TEMPORARY BIPOLAR BALLOON PACING CATHETER FLOW DIRECTED


- intravenous placement without fluoroscopic monitoring
- allows for placement at bedside
- for pre- and postoperative and emergency stimulation
- interim support during pacemaker exchange
- intracardial monitoring

FEATURES

- unique design: particularly soft, flexible tip with solid catheter body for excellent torqueability
- antithrombogenic polyurethane catheter body
- available with shrouded (touch protected) pins
- unshrouded connectors are only for use with unshrouded 2mm receptacles at external pacemakers
- radiopaque
- 10 cm depth markings facilitate precise and safe positioning



shrouded or unshrouded connectors



flexible soft tip



depth markings



flow directed bipolar pacing catheter

TEMPORARY BIPOLAR PACING CATHETER BIPOLAR ELECTRODE FOR TEMPORARY INTRACARDIAC PACING

FIELDS OF APPLICATION:

- Cardiology
- Radiology
- Anaesthesia
- Emergency

CONNECTOR CAPS:

Supplied separately with the catheter for touch protected connection

PACKAGING:

double packed, sterile

MATERIAL:

Catheter tube: PUR, radiopaque

Electrodes: stainless steel

Balloon: Latex

DIMENSIONS:

Usable length: 110 cm

Depth markings: every 10 cm

Length of flexible distal part: 70 mm

Number of electrodes: 2

Electrode distance: 10 mm

TEMPORARY BIPOLAR PACING CATHETER, TORQUE DIRECTED

Pacing Catheter, FEP Introducer Cannula

| Ref.No. Standard Tip | Ref.No. Soft Tip | French Size | Electrode Curve | Rcmd. Introducer Size | Max. Balloon Diameter | Packing Unit |
|-------------------------|---------------------|----------------|--------------------|-----------------------------|-----------------------------|-----------------|
| 024.TPL4 | - | 4 F / 1.3 mm | Cournand | 5 F | - | 5 each |
| 024.TPL5 | 025.TPL5 | 5 F / 1.7 mm | Cournand | 5 F | - | 5 each |
| 024.TPL6 | 025.TPL6 | 6 F / 2.0 mm | Cournand | 6 F | - | 5 each |

Pacing Catheter, without Introducer Cannula

| | | | | | | |
|-----------|-----------|--------------|----------|-----|---|--------|
| 024.TPL70 | 025.TPL70 | 7 F / 2.3 mm | Cournand | 7 F | - | 5 each |
|-----------|-----------|--------------|----------|-----|---|--------|

Available on request:

- Pacing catheter with J-tip or straight tip
- Pacing catheter sizes 5F and 6F with Peel-away cannula

TEMPORARY BIPOLAR PACING CATHETER

BIPOLAR ELECTRODE FOR TEMPORARY INTRACARDIAC PACING

BIPOLAR PACING CATHETER SETS, TORQUE DIRECTED

Pacing Catheter with Introducer Set:
Dilator 12 cm / Sheath 6 cm, Seldinger Cannula 18G, Guidewire (double ended) .035" x 45 cm

| Ref.No. Standard Tip | Ref.No. Soft Tip | French Size | Electrode Curve | Rcmd. Introducer Size | Max. Balloon Diameter | Packing Unit |
|-------------------------|---------------------|----------------|--------------------|-----------------------------|-----------------------------|-----------------|
| 024.TPLS14 | - | 4 F / 1.3 mm | Cournand | 5 F | - | 5 each |
| 024.TPLS15 | - | 5 F / 1.7 mm | Cournand | 5 F | - | 5 each |
| 024.TPLS16 | - | 6 F / 2.0 mm | Cournand | 6 F | - | 5 each |
| 024.TPLS17 | - | 7 F / 2.3 mm | Cournand | 7 F | - | 5 each |

Pacing Catheter with Introducer Set:
Dilator 12 cm / Sheath 6 cm, Touhy-Borst Valve, Seldinger Cannula 18G,
Guidewire (double ended) .035" x 45 cm, Sterility Sleeve

| Ref.No. Standard Tip | Ref.No. Soft Tip | French Size | Electrode Curve | Rcmd. Introducer Size | Max. Balloon Diameter | Packing Unit |
|-------------------------|---------------------|----------------|--------------------|-----------------------------|-----------------------------|-----------------|
| 024.TPLS24 | - | 4 F / 1.3 mm | Cournand | 5 F | - | 5 each |
| 024.TPLS25 | 025.TPLS25 | 5 F / 1.7 mm | Cournand | 5 F | - | 5 each |
| 024.TPLS26 | 025.TPLS26 | 6 F / 2.0 mm | Cournand | 6 F | - | 5 each |
| 024.TPLS27 | 025.TPLS27 | 7 F / 2.3 mm | Cournand | 7 F | - | 5 each |

BIPOLAR BALLOON PACING CATHETER, FLOW DIRECTED

Balloon Pacing Catheter, Syringe

| Ref.No. Standard Tip | French Size | Electrode Curve | Rcmd. Introducer Size | Max. Balloon Diameter | Packing Unit |
|-------------------------|----------------|--------------------|-----------------------------|-----------------------------|-----------------|
| 024.BP51 | 5 F / 1.7 mm | Straight | 6 F | 9 mm | 1 each |

BIPOLAR BALLOON PACING CATHETER SETS, FLOW DIRECTED

Balloon Pacing Catheter with Introducer Set:
Dilator 12 cm / Sheath 6 cm, Touhy-Borst Valve, Seldinger Cannula 18G,
Guidewire (double ended) .035" x 45 cm, Sterility Sleeve

| Ref.No. Standard Tip | French Size | Electrode Curve | Rcmd. Introducer Size | Max. Balloon Diameter | Packing Unit |
|-------------------------|----------------|--------------------|-----------------------------|-----------------------------|-----------------|
| 024.BPS51 | 5 F / 1.7 mm | Straight | 6 F | 9 mm | 5 each |